Exacto®

FR

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REF 0075N

BIOSYNEX
EASY DIAGNOSTICS FOR LIFE

SOMMAIRE

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
BIOSYNEX SA
22 Boulevard Sébastien Brant
67400 Illkirch-Graffenstaden
FRANCE

Date de dernière modification : 01/2020

 AViTA Corporation
9F, NO.78, Sec.1, Kwang-Fu
Rd., san-Chung District 24158,
New Taipei City,
Taiwan, R.O.C.

client@biosynex.com
www.biosynex.com

 MDSS :
Medical Device Safety Service GmbH (MDSS)
Schiffgraben 41,
D-30175 Hannover,
Germany

CE 0598

IFU_0075N_ML_V07202001R02

Exacto®

EN

Blood pressure monitor

REF 0075N

BIOSYNEX
EASY DIAGNOSTICS FOR LIFE

SUMMARY

Introduction	23
Blood Pressure Monitor Features	23
Important Information before use	24 - 25
LCD screen description	25
Battery installation	26
Setting the date and time	27
Setting the alarm clock	28
Placement of the pressure sleeve	29 - 30
Measurement of pulse rate and blood pressure	31 - 32
Irregular Heartbeat Detector	33
World Health Organization (WHO)	33
Memory Function	34
Troubleshooting	35
Error Code	36
Care and Maintenance	37
Standards applied	37
Technical Specifications	38
EMC Tables	38 - 39
Symbols Index	40

INTRODUCTION

To get the maximum benefit from your blood pressure monitor, we recommend that you consult your doctor or qualified health care professional first.

Thank you for buying a blood pressure monitor. When properly used and carefully maintained, your blood pressure monitor will provide you with years of reliable readings.

The measurement method used by your automatic blood pressure monitor is called the oscillometric method. The blood pressure monitor detects your blood's movement in the artery and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the blood pressure monitor easier to use.


The pre-shaped cuff is easy to put on and comfortable to wear.

Frequently measuring blood pressure will provide you and your doctor with an accurate indication of changes in your actual blood pressure.

Those with severe hypertension, severe arteriosclerosis or severe diabetes must consult their doctor regarding blood pressure monitoring in the arm.

Blood pressure measurements with this device are equivalent to measurements obtained by a trained observer using the cuff/stethoscope auscultation method within the limits prescribed by the National Standard for Electronic or Automated Sphygmomanometers.

BLOOD PRESSURE MONITOR FEATURES

IMPORTANT INFORMATION BEFORE USE

When you use this product, be sure to follow all instructions below. Any action contrary to these instructions may affect the accuracy of the device.

1. Do not modify this equipment without the authorization of the manufacturer.
2. Blood pressure measurements should only be interpreted by a qualified medical professional who knows your medical history.

By using this device regularly and recording your measurements, you can keep your doctor informed of changes in your blood pressure

3. Perform your measurements in a quiet place.

You should be seated in a relaxed position.

4. Avoid smoking, eating, taking medication, consuming alcohol, and any physical activity 30 minutes before taking your blood pressure.

If you are showing signs of stress, avoid taking your blood pressure until the feeling has faded.

5. Rest for 15 minutes before taking your blood pressure.
6. Take off tight clothing or jewellery that may interfere with the placement of the cuff.
7. Keep the blood pressure monitor stable while measuring to get an accurate reading.

Keep still and do not speak during the measurement.


8. Record your daily blood pressure and pulse readings.
9. Always take your blood pressure at the same time each day or according to your doctor's advice to get an accurate indication of your actual changes in blood pressure.
10. Wait at least 15 minutes between measurements to allow blood vessels to return to normal.
11. In the event that your pulse is very weak or irregular, errors may occur that prevent correct measurement. If any abnormal changes are observed, consult your doctor or health care professional.

12. This device is intended for adult use. When taking a measurement, you can stop the cuff inflation or deflation at any time by pressing the POWER button.

13. If you need help during the setup, use or care, please contact the manufacturer.

Using a blood pressure monitor cannot in any case exempt you from consulting a doctor. Do not change your treatment without medical advice.

LCD SCREEN DESCRIPTION


BATTERY INSTALLATION

Low battery warning: 

To change the batteries:

1. Open the battery cover.
2. Making sure the battery polarity is correct, install four 1.5 V AA batteries.
3. Close the battery cover.


NOTE:

1. Please properly dispose of the batteries away from small children and heat.
2. It is recommended to remove the batteries if the unit will not be used for an extended period of time.
3. Batteries must be disposed of in accordance with local environmental and institutional policies.

Setting the date and Time

It is necessary to set the date and time for the unit every time batteries are initially installed or replaced.


1. While in power off mode, press the "Time Set" button to enter Date and Time setting procedure and the Year value will begin to flash.
2. Press the "MEMORY Recall" button to advance the display to the desired year, press "POWER" button to confirm the year.
3. Next, the Month will blink. Repeat steps 2 to set the month and date, then hours, then minutes.
4. After setting the minutes, the unit will automatically exit out of the date/time setting mode and briefly show the word OFF before shutting down.

SETTING THE ALARM CLOCK


This monitor features an alarm clock setting, allowing the user to set reminders for up to three daily alarms. The alarm clock function serves as a reminder for those users who may need to track measurements at consistent, specific times of the day.

1. Press the " Alarm set" button to enter Alarm Clock Setting procedure.

2. Continue to press the "MODE" button until you have selected which alarm to set (alarm 1, 2 or 3).


3. With the chosen alarm screen displayed, the hour will begin to flash.


4. Press the "MEMORY Recall" button to advance to the desired "hour".

Press the " POWER" button to confirm.

5. Press the "MEMORY Recall" button to advance to the desired "minutes". Press the " POWER" button to confirm.

6. Press the "MEMORY Recall" button to toggle between ON and OFF. Press the " POWER" button to confirm your selection.

7. If setting more than one alarm reminder, repeat all steps for the 2nd and 3rd alarm clock.

8. Press the POWER button to turn off the unit. The alarm will sound at the appropriate set time even though the unit is powered OFF.

PLACEMENT OF THE PRESSURE SLEEVE

It is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading. If for any reason you are unable to or should not use your left arm, please modify the instructions for cuff application to your right arm. Your physician can help you identify which arm is best for you to take measurements from.

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.




2. Be seated at a table or desk with your feet flat on the floor.

3. The cuff should not be plugged into the monitor until after the cuff is applied to your arm.


4. Position the cuff on a solid surface with the tubing facing up and away from you. The metal ring/bar on the cuff should be to the left of the tubing.

5. Open the cuff by pulling or rolling the bottom of the cuff to the right. This should open the cuff without fully unrolling it, creating a cylinder. Do not fully unwrap or unroll the cuff.


6. Insert your left arm into the created cuff cylinder. Position the (♣) mark over the main artery on the inside of your arm.

7. The bottom edge of the cuff should be positioned approximately one inch above the elbow joint.

8. Reaching underneath your left arm with your right hand, pull the end of the cuff towards your body to tighten the cuff. Wrap and secure the cuff, making sure that the (Ψ) mark remains in place as shown.

9. The cuff should fit comfortably, yet snugly around your arm. You should be able to insert one finger easily between your arm and the cuff.

10. The cuff size is suitable for use when the vertical 'index' mark is within the horizontal OK range. A different size cuff is needed if the 'index' mark is outside the OK range.

If you are not comfortable with applying your cuff, please seek the assistance of another member of your household or work with your physician to practice the cuff application. Incorrectly applied cuffs may result in inaccurate readings.

MEASUREMENT OF PULSE RATE AND BLOOD PRESSURE

Please read the preceding portions of this manual prior to taking your first reading.

1. Position the monitor on a flat, stable surface with the digital display in view.

2. Insert the cuff tubing connector into the port on the left side of your monitor.

3. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart. Relax your left hand.




4. Press the POWER button. This will turn the power on.

5. After the self-test, the values for the last reading will appear on the display.


6. The blood pressure monitor start to measure.

7. The cuff will automatically begin to inflate, with the display showing the increasing pressure in the cuff. As the pressure increases, an arrow pointing up will appear on the display.

8. When the inflation has reached optimum level, the display will begin to show the decreasing pressure; the screen will display an arrow pointing down while you feel the pressure of the cuff decrease.

9. To detect the heartbeat, the heartbeat symbol will appear and continuous flashes on the LCD display.

10. Your blood pressure measurement and pulse will display simultaneously on the screen.

11. The Hypertension Indicator will indicate your reading range on the display separately.

12. Your reading will be automatically stored in «MEMORY Recall» Button.

13. Select POWER to turn the unit off and conserve energy and battery life. The unit will automatically shut-off approximately 2 minutes.

IHD: IRREGULAR HEARTBEAT DETECTOR

1. Your digital blood pressure monitor has an irregular heartbeat detector (IHD). An irregular heartbeat may influence the measurement. If the blood pressure monitor detects an irregular heartbeat during the measurement, the «» icon will appear on the screen with the measurement values. You may take another measurement to ensure that the values are not being influenced by movement during the reading or by an irregular heartbeat.


2. If this icon appears frequently, we advise you to consult with your doctor.

3. This device does not replace a cardiac exam, but serves to detect heartbeat irregularities at an early stage.

WORLD HEALTH ORGANIZATION (WHO)

This device is equipped with our unique hypertension indicator. The World Health Organization has established globally recognized standards to evaluate high or low blood pressure readings. The graphic below should only be considered a guideline; always consult your doctor or health care professional to interpret your individual results.


MEMORY FUNCTION

Recall measurements in the Memory

You can recall up to 60 total measurements in the memory, plus an average of the 3 last measurements saved in the memory to share them with your doctor or a qualified health care professional.

1. Press and release the «MEMORY» button. The device will then show the average of the last 3 saved measurements.
2. Continue to push the «MEMORY» button to see the previously saved measurements. The measurements will appear on the screen from the most recent to the oldest; the number of memory will appear on the lower left corner.

3. All the results of a given measurement will show, including the measurement results, the heart rate, the hypertension indicator and the timestamp.
4. When the number of readings surpasses 60, the oldest data will be replaced by the new recording.
5. Press the Power button to turn off the blood pressure monitor at any time while reading saved measurements.

Erasing Measurements from the Memory:

From the turned off screen, hold the «MEMORY» button down until the screen shows CLr. This indicates that all measurements have been erased.

TROUBLESHOOTING

Issue	Likely cause	Recommended action
Nothing appears on the screen even when the device is powered on.	The batteries are dead.	Replace all batteries with new batteries.
	The batteries are not properly aligned with the tracks.	Replace the batteries in the proper position.
The low battery icon appears.	The batteries are dead.	Replace all batteries with new batteries.
	In very cold conditions, the batteries have lower electrical charges.	Warm up the batteries or use the device in a warmer environment.
The device's operating time is inconsistent.	Different brands of batteries have different lifetimes.	Use alkaline batteries and replace all the batteries at the same time with the same brand of batteries.
No reading after the measurement.	The batteries are dead.	Replace all batteries with new batteries.
Questionable blood pressure results.	The cuff may have been poorly positioned.	Adjust the cuff to measure.
	Blood pressure varies naturally throughout the day.	Rest for a moment, relax and measure again.
Questionable heart rate results.	Body movement during use of the device.	Avoid moving while taking the reading.
	Measurements taken shortly after exercise or time outside.	Do not take measurements after exercise or coming back in from outside.
The device turns off automatically.	Device design.	Press the power button again, then measure again.
During the measurement, the cuff begins to inflate again.	This could be a normal action if the user's blood pressure is higher than the initial pressure value, the device automatically pumps to a higher pressure level of 40mmHg each time.	Relax, and try to take a measurement again.
	The cuff is not secured.	Make sure that the arm cuff is securely fastened and retake the measurement.

ERROR CODE

Error code	Meaning	Corrective action
Err 0	No pulse or insufficient pulse detection	Remove heavy clothing and retry.
Err 1	Leak in the pressure cuff/ inflation too weak	The arm cuff is not securely fastened. Put it on again and retake the measurement.
Err 2	Pressure failure	Rest for a moment, relax and retry.
Err 3	Deflation failure	The arm cuff is not securely fastened. Put it on again and retake the measurement.
Err	Memory error	Take out the batteries to restart the device, then retake the measurement.
	Low batteries	Replace all batteries with new batteries.

CARE AND MAINTENANCE

1. Clean the device and arm cuff carefully only with a slightly moistened cloth.
2. Do not immerse the device in water. It may cause damage if water enters.
3. Do not use any gas, strong detergent, or solvents to clean the device, including the cuff.
4. Disconnect the cuff and tubing from the monitor prior to storing.
5. Do not use any liquids on the monitor or cuff.
6. Use a soft, dry cloth to clean your monitor.
7. Do not store in direct sunlight, dust or excessive humidity.
8. Avoid extreme temperatures.
9. Do not disassemble the monitor or cuff.
10. Remove batteries if the monitor will not be used for an extensive period of time.

STANDARDS APPLIED

This product is in accordance with the measures in directive CE MDD (93/42/CEE).
The following standards apply to the design and/or the manufacturing of the products:

EN 1060-3

Non-invasive sphygmomanometers – Part 3: Additional requirements for the electromechanical blood pressure measuring system

EN 1060-4

Non-invasive sphygmomanometers – Part 4: Test procedures that aim to determine the overall system accuracy of automated non-invasive sphygmomanometers

ISO 14971

Medical Devices – Application of risk management to medical devices. Classification according to IEC/EN 60601-1 paragraph 5:

- Internally-powered equipment
- IPX0
- Equipment not suitable for use in the presence of an anaesthetic mixture flammable with air or oxygen or nitrous oxide.
- Continuous operation

TECHNICAL SPECIFICATIONS

• Measuring range:	Blood Pressure : 30 à 280 mmHg Pulse rate : 40 à 199 pulse / min
• Calibration Accuracy;	Blood Pressure : ± 3 mmHg Pulse rate: $\pm 4\%$ of reading
• Operating Environment	10°C à 40°C with relative humidity up to 85% (non condensing)
• Storage / Transport Environment	-20 °C à + 50 °C with relative humidity up to 85% (non condensing)
• Power source:	4 x AA 1,5 V (LR06)
• Weight:	259g (without batteries)
• Dimensions :	128 mm x 150 mm x 59,2 mm (l x H x P)
• Cuff Circonférence (M-L size) :	22 - 42 cm

EMC TABLES

Manufacturer's declaration and directive – electromagnetic emissions

The 0075N arm blood pressure monitor is intended for use in the electromagnetic environment specified below. The client or user of the 0075N arm blood pressure monitor must be sure to use it in such an environment.

Emissions test	Compliance	Electromagnetic environment guidelines
RF Emissions CISPR 11	Group 1	The 0075N arm blood pressure monitor uses RF energy only for its internal operation. Therefore, its RF emissions are very low and are not likely to cause interference in electronic equipment nearby.
RF Emissions CISPR 11	Class B	
Harmonic Emissions CEI 61000-3-2	Not applicable	
Voltage variations / Flicker emissions CEI 61000-3-3	Not applicable	

EMC TABLES

Manufacturer's declaration and directive – electromagnetic immunity

The 0075N arm blood pressure monitor is intended for use in the electromagnetic environment specified below. The client or user of the 0075N arm blood pressure monitor must be sure to use it in such an environment.

Immunity Test	CEI 60601 Test level	Compliance level	Guidelines on the electromagnetic environment
Electrostatic Discharge (DES) CEI 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	The floors must be made of wood, cement or ceramic tiles. If the floors are covered in synthetic material, relative humidity must be below 30%.
Power Supply Frequency (50/60 Hz) Magnetic Field CEI 61000-4-8	3 A/m	3 A/m	The magnetic field of the supply frequency must be at a level appropriate to the place of use.

EMC TABLES

Manufacturer's declaration and directive – electromagnetic immunity

The 0075N arm blood pressure monitor is intended for use in the electromagnetic environment specified below. The client or user of the 0075N arm blood pressure monitor must be sure to use it in such an environment.

Immunity Test	CEI 60601 Test level	Compliance level	Electromagnetic environment guidelines
Radiated RF CEI 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Portable and mobile RF communication equipment must not be used closer to any part of the 0075N arm blood pressure monitor, including the cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $D = 1.2 \sqrt{P}$ $D = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $D = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). The field strengths of fixed RF transmitters, as determined by an electromagnetic site survey, must be less than the compliance level in each frequency range. Interference may occur near equipment marked with the following symbol: </p>

NOTE 1 At 80 MHz and 800 MHz, the highest frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.

A. The field strengths of stationary transmitters, such as base stations for radio (cellular/wireless) and land mobile radios, amateur radios, AM and FM broadcasting and TV broadcasting may not in theory be predicted accurately. To evaluate the electromagnetic environment due to fixed RF transmitters, consideration should be given to an electromagnetic survey of the site. If the measured field strength value at the location of the 0075N arm blood pressure monitor exceeds the applicable RF compliance level above, the 0075N arm blood pressure monitor should be checked to verify normal operation. If abnormal performance is observed, additional measures may be required, such as reorienting or reinstalling the 0075N arm blood pressure monitor.

B. Depending on the frequency range from 150 kHz to 80 MHz, the field strengths must be less than 3 V/m.

EMC TABLES

Recommended separation distance between portable and mobile RF communication equipment and the 0075N arm blood pressure monitor.

The 0075N arm blood pressure monitor is intended for use in an electromagnetic environment where the RF perturbations are controlled. The customer or user of the 0075N arm blood pressure monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications (Transmitters) and the 0075N arm blood pressure monitor as recommended below, depending on the maximum output power of the communication equipment.

Maximum rated output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters whose maximum output power is not specified above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the Maximum output power of the transmitter in watts (W) according to the manufacturer of the transmitter.

NOTE 1 At 80 MHz and 800 MHz, separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.

SYMBOLS INDEX

	Read the instruction manual before use
	Type BF applied part
	Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations
	Manufacturer
	Storage conditions
REF	Catalog number
	Keep dry
LOT	Lot number
CE 0598	This device complies with the requirements of the Medical Device Directive 93/42/EEC.
	Common signage TRIMAN
	Production year
EC REP	European Authorized Representative

WARRANTY

This product is warranted for a two-year period from the purchase date. However, this warranty does not apply to damages, in particular, resulting from non-compliance with this guide, incorrect manipulations, shocks or inappropriate installation of the battery. It does not apply either to the battery or the package.

Imported by :

BIOSYNEX SA
22 Boulevard Sébastien Brant
67400 Illkirch-Graffenstaden
FRANCE

 AViTA Corporation
9F, NO.78, Sec.1, Kwang-Fu
Rd., san-Chung District 24158,
New Taipei City,
Taiwan, R.O.C.

client@biosynex.com
www.biosynex.com

Date of Last Modification : 01/2020

EC **REP**

MDSS:
Medical Device Safety Service GmbH (MDSS)
Schiffgraben 41,
D-30175 Hannover,
Germany

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AU S.A.V. BIOSYNEX L'ADRESSE INDIQUEE CI-DESSOUS